

ERVIN DANESH AZMA



ERVIN DANESH AZMA

# اروین دانشن آزما

Ervin Danesh Azma

✓ خدمات خوب یک اتفاق نیست





اروین دانش آزما

Ervin Danesh Azma

✓ خدمات خوب یک اتفاق نیست



# MYTHIC™ 5VET PRO

Fully automated 5DIFF hematology analyzer for veterinary use.





اروین دانش آزما  
Ervin Danesh Azma

✓ خدمات خوب یک اتفاق نیست

## Technical specifications

|                           |  |
|---------------------------|--|
| Throughput                | 60 samples per hour  |
| Sample volume             | 20 µL (Whole blood and pre-diluted mode)   |
| Reagents                  | Diluent, Optical Fluid, Lysing Reagent, Cleaner reagents recognized with RFID code   |
| Technology                | <ul style="list-style-type: none"> <li>• Laser Light Multi-Dimensional Cell</li> <li>• Classification (WBC DIFF)</li> <li>• Impedance (WBC, RBC, PLT)</li> <li>• Spectrophotometry (HGB)</li> </ul>        |
| 34 parameters             | WBC, LYM% and #, MON% and #, NEU% and #, EOS% and #, BASO% and #, RBC, HGB, HCT, MCV, MCH, MCHC, RDW_CV, RDW_SD, PLT, MPV, PDW, PCT, P_LCC, P_LCR, RET%, RET_ABS, IRF, NRBC% and #, LIC% and #, ALY% and # |
| Histograms & Scattergrams | RBC histogram, PLT histogram, Scattergrams for WBC differentiation, Scattergrams for RET   |
| Dimensions                | 270 x 508 x 412 mm (W x D x H)   |
| Weight                    | 26 kg  |
| Power                     | <ul style="list-style-type: none"> <li>• 100 V to 240 V</li> <li>• 50 to 60 Hz</li> <li>• Less than 250 VA</li> </ul>  |
| Printers                  | <ul style="list-style-type: none"> <li>• Built-in thermal printer</li> <li>• External printer USB port</li> </ul>  |
| User's interface          | <ul style="list-style-type: none"> <li>• 10.4" touchscreen</li> <li>• info about keyboard to be deleted</li> </ul>   |
| Output                    | <ul style="list-style-type: none"> <li>• RS232</li> <li>• Ethernet (UDP/IP, TCP/IP protocol)</li> <li>• 2 USB port</li> <li>• Centronics port</li> <li>• Barcode scanner port</li> </ul>                   |

## Software specifications

|                   |  |
|-------------------|--|
| Graphic interface | Intuitive menu with big icons  |
| Available species | Dog, Cat, Horse, Pig, Cow, Buffalo, Rabbit, Monkey, Rat, Mouse, Sheep, Goat, Camel, 3 User-defined   |
| Multilanguage     | English, Chinese, Spanish, French, Russian, Portuguese, Czech, Turkish, Polish, Italian  |
| Calibration       | Automatic or manual  |
| Quality control   | <ul style="list-style-type: none"> <li>• 6 batches: 100 run/batch</li> <li>• Table and Levey Jennings graphs</li> </ul>  |
| Archives          | <ul style="list-style-type: none"> <li>• Internal memory: 200 000 patients</li> </ul>  |
| Traceability      | <ul style="list-style-type: none"> <li>• Patient identification: name, patient ID, sample ID</li> <li>• Daily log</li> <li>• Reagents monitoring: lot number, date, volume</li> <li>• Service log</li> </ul> |
| Flagging          | <ul style="list-style-type: none"> <li>• Pathologic flags</li> <li>• Lab limits</li> <li>• Reagent alerts</li> <li>• Instrument alerts</li> </ul>  |
| Reagents          | Mythic 5 VET PRO Diluent<br>Mythic 5 VET PRO Optical Fluid<br>Mythic 5 VET PRO Lysing Reagent<br>Mythic 5 VET PRO Cleaner<br>MYTHIC RETIC PRO  |



Professional Veterinary use Hematology system for the complex 5DIFF blood count with Reticulocytes determination

- 1 Orphée reserves own rights to technical and software modification. In case of printer choice please contact your local Orphée's representative. Number of tests per pack depends on the analyzed species and daily workflow. Mass storage depends on capacity of USB mass storage device.

# اروین دانش آزما

Ervin Danesh Azma



021\_62999800

09100145015\_18



CO.ERVINDANESHAZMA



WWW.ERVINDANESHAZMA.COM